

## **Neck Collar**

For people with limited head control





Children, adolescents, and adults with cerebral palsy (CP) or with other types of brain damages, syndromes or disabilities where the muscles are too weak to support the weight of the head, often have challenges holding the head in an upright position.

Insufficient control of the head due to limited neck muscle control and strength affects positioning in sitting, standing and walking. Inappropriate positioning may influence the person's ability to communicate and interact with others, as well as having negative effects on respiration and digestion.

For this reason, appropriate positioning of the head has become an accepted goal in an upright sitting or standing position.

Made for Movement offers a Neck Collar that positions the persons head in a natural head position in both sitting, standing and walking situations.

The Neck Collar provides necessary support, without the user being deprived of full head control.

## **Advantages of the Neck Collar**

- Easy to put on and take off
- Provides necessary support of the head in a natural position
- · Allows movement of the head
- · Can be used during transfer
- Enhances eye to eye contact

model	S	М	L
Circumference (+/- 5 cm)	30 cm	40 cm	50 cm

The Neck collar comes in three sizes and is easily adjustable with a Velcro strap. This means that the Neck Collar will suit most users, both children and adults.

They are designed to be comfortable to wear. The cover is removable and can be washed and disinfected.



## The uses are many, and the Made for Movement's Neck Collar provides good stability and support without restriction.









In transfer situations

In wheelchair/seated

In gait trainer

In standing/ dynamic stander

